

CEL-Seq2 Protocol in the C1™

This document is a merging of Supplementary file 1 and Supplementary file 2 from the publication T. Hashimshony, T. N. Senderovich, G. Avital, A. Klochendler, Y. de Leeuw, L. Anavy, D. Gennery, S. Li, K. J. Livak, O. Rozenblatt-Rosen, Y. Dor, A. Regev, and I. Yanai. CEL-Seq2: sensitive highly-multiplexed single-cell RNA-Seq. *Genome Biol* **17**:77 (2016).

Reagents:

Ultra pure RNase free water

Ethanol

Bioanalyzer kits - Agilent RNA pico kit (5067-1513), high sensitivity DNA kit (5067-4626)

Qubit reagents: dsDNA HS Assay – Invitrogen Q32851 or Q32854

For RNA amplification:

ERCC RNA spike-in mix – Thermo Fisher Scientific 4456740

C1 Open App IFC – Fluidigm 100-8133 (5–10 µm), 100-8134 (10–17 µm), 100-8135 (17–25 µm)

C1 Open App Reagent Kit – Fluidigm 100-8920

NP-40, 10% – Thermo Fisher Scientific 28324

5× Lysis buffer: 2.5% NP-40, 250 mM Tris-HCl, pH 8.4, 5 mM EDTA

DTT (dithiothreitol) – Teknova D9750

dNTP Mix, 25 mM each – Thermo Fisher Scientific R1121

rNTP Mix, 25 mM each – New England BioLabs N0466S

SuperScript II – Thermo Fisher Scientific 18064014 (includes 5× First Strand Buffer)

RNaseOUT – Thermo Fisher Scientific 10777019

Second Strand Buffer – Thermo Fisher Scientific 10812014

DNA Polymerase I (E. coli) – Thermo Fisher Scientific 18010025

E. coli DNA ligase – Thermo Fisher Scientific 18052019

RNaseH (E. coli) – Thermo Fisher Scientific 18021071

MEGAscript T7 Transcription Kit – Thermo Fisher Scientific AM1334

ExoSAP-IT for PCR Product Clean-Up – Affymetrix 78200

Fragmentation buffer: 200 mM Tris-acetate, pH 8.1, 500 mM KOAc, 150 mM MgOAc

Fragmentation stop buffer: 0.5 M EDTA pH 8

AMPure XP beads – Beckman Coulter A63880

RNAClean XP beads - Beckman Coulter A63987

Bead binding buffer - 20% PEG8000, 2.5 M NaCl

For Library preparation:

SuperScript II – Thermo Fisher Scientific 18064014

RNaseOUT – Thermo Fisher Scientific 10777019

AMPure XP beads – Beckman Coulter A63880

Phusion® High-Fidelity PCR Master Mix with HF Buffer – New England BioLabs M0531

randomhexRT primer – GCCTGGCACCCGAGAATTCCANNNNNN

RNA PCR primers (sequences available from Illumina)

RP1 – AATGATACGGCGACCACCGAGATCTACACGTTTCAGAGTTCTACAGTCCG*A, *=phosphorothioate

RPi1 – CAAGCAGAAGACGGCATAACGAGATCGTGATGTGACTGGAGTTCCTTGGCACCCGAGAATTCC*A

Primers re-suspended at 100 µM, PCR primers used at 10 µM.

Equipment:

Thermocycler with lid with adjustable temperature (one that can also fit 0.5 ml PCR tubes is convenient)

Oven (optional)

Magnetic stand (for 0.5 ml tubes)

Qubit® Fluorometer - invitrogen

Bioanalyzer – Agilent

Primers:

CEL-Seq2 primer design: The RT primer was designed with an anchored polyT, a 6 bp unique barcode, a 6 bp UMI (unique molecular identifier), the 5' Illumina adapter (as used in the Illumina small RNA kit) and a T7 promoter. The barcodes were designed such that each pair is different by at least two nucleotides, so that a single sequencing error will not produce the wrong barcode. Primers are desalted at the lowest possible scale, stock solution 1 µg/µl, working concentration 25 ng/µL.

1s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGACTCTTTTTTTTTTTTTTTTTTTTTT
2s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGCTAGTTTTTTTTTTTTTTTTTTTTT
3s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGCTCATTTTTTTTTTTTTTTTTTTTTT
4s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGCTTCTTTTTTTTTTTTTTTTTTTTTT
5s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCATGAGTTTTTTTTTTTTTTTTTTTTT
6s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCATGCATTTTTTTTTTTTTTTTTTTTTT
7s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCATGTCTTTTTTTTTTTTTTTTTTTTTT
8s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCACTAGTTTTTTTTTTTTTTTTTTTTT
9s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCAGATCTTTTTTTTTTTTTTTTTTTTTT
10s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCACAGTTTTTTTTTTTTTTTTTTTTT
11s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGGATCTTTTTTTTTTTTTTTTTTTTTT
12s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGTGATTTTTTTTTTTTTTTTTTTTTT
13s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNAGTGTCTTTTTTTTTTTTTTTTTTTTTT
14s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTAGTTTTTTTTTTTTTTTTTTTTT
15s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTGAGTTTTTTTTTTTTTTTTTTTTT
16s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTGCATTTTTTTTTTTTTTTTTTTTTT
17s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCGAAGTTTTTTTTTTTTTTTTTTTTT
18s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCGACATTTTTTTTTTTTTTTTTTTTTT
19s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCGATCTTTTTTTTTTTTTTTTTTTTTT
20s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGTACAGTTTTTTTTTTTTTTTTTTTTT
21s	GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGTACCATTTTTTTTTTTTTTTTTTTTTT

67s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCACCATTTTTTTTTTTTTTTTTTTTTTV
68s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCACTCTTTTTTTTTTTTTTTTTTTTTTV
69s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCCTCATTTTTTTTTTTTTTTTTTTTTTV
70s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTCCTTCTTTTTTTTTTTTTTTTTTTTTTV
71s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTGTCTTTTTTTTTTTTTTTTTTTTTTV
72s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGTCTTCTTTTTTTTTTTTTTTTTTTTTTV
73s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGTTGAGTTTTTTTTTTTTTTTTTTTTTV
74s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGTTGTCTTTTTTTTTTTTTTTTTTTTTTV
75s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGTGAAGTTTTTTTTTTTTTTTTTTTTTV
76s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNACAGACTTTTTTTTTTTTTTTTTTTTTTV
77s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNACAGGATTTTTTTTTTTTTTTTTTTTTTV
78s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNACCAACTTTTTTTTTTTTTTTTTTTTTTV
79s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNACCAGATTTTTTTTTTTTTTTTTTTTTTV
80s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNACTCACTTTTTTTTTTTTTTTTTTTTTTV
81s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTCAACTTTTTTTTTTTTTTTTTTTTTTV
82s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTTCACTTTTTTTTTTTTTTTTTTTTTTV
83s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTTCTGTTTTTTTTTTTTTTTTTTTTTV
84s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNCTGTTGTTTTTTTTTTTTTTTTTTTTTV
85s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTGAGTGTTTTTTTTTTTTTTTTTTTTTTV
86s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTGAGGATTTTTTTTTTTTTTTTTTTTTTV
87s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTGTCTGTTTTTTTTTTTTTTTTTTTTTV
88s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTGGTTGTTTTTTTTTTTTTTTTTTTTTV
89s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNTGGTGATTTTTTTTTTTTTTTTTTTTTTV
90s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGAAGACTTTTTTTTTTTTTTTTTTTTTTV
91s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGAAGTGTTTTTTTTTTTTTTTTTTTTTTV
92s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGAAGGATTTTTTTTTTTTTTTTTTTTTTV
93s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGACAACTTTTTTTTTTTTTTTTTTTTTV
94s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGACAGATTTTTTTTTTTTTTTTTTTTTTV
95s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGAGTTGTTTTTTTTTTTTTTTTTTTTTV
96s GCCGGTAATACGACTCACTATAGGGAGTTCTACAGTCCGACGATCNNNNNNGAGTGATTTTTTTTTTTTTTTTTTTTTTV

C1 Mixes

Lysis Mix (1000 µL)

280 µL 5× Lysis Buffer
28 µL 1 M DTT
4 µL 40 units/µL RNaseOUT
8 µL 1:10,000 dilution of ERCC RNA Spike-In Mix
70 µL C1 Loading Reagent
610 µL H₂O
7 µL Lysis Mix to Inlet #3

Diluted CEL-Seq2 Primers

Dispense 8 µL Lysis Mix to each of 96 wells
Add 2 µL 25 ng/µL CEL-Seq2 Primer to each well
5 µL Diluted CEL-Seq2 Primer to each of 96 Harvest Outlets

RT Mix (28.6 µL)

20 µL 5× First Strand Buffer
2.6 µL 40 units/µL RNaseOUT
2 µL 25 mM each dNTP
2.6 µL 200 U/µL SuperScript II Reverse Transcriptase
1.4 µL C1 Loading Reagent
7 µL RT Mix to Inlet #4

2nd Strand Mix (100 µL)

24.7 µL 5× Second Strand Buffer
1 µL 25 mM each dNTP
3.3 µL 10 units/µL E.coli DNA polymerase
1 µL 10 units/µL E. coli DNA ligase
1 µL 2 units/µL RNase H
5 µL C1 Loading Reagent
64 µL H₂O
24 µL 2nd Strand Mix to Inlet #7

IVT Mix (29.6 µL)

6.6 µL T7 10× Reaction Buffer (MEGAscript kit)
14.9 µL 25 mM each rNTP
6.6 µL T7 Enzyme Mix (MEGAscript kit)
1.5 µL C1 Loading Reagent
24 µL IVT Mix to Inlet #8

C1 Steps

1. 96 Diluted CEL-Seq2 Primers to chamber 1 from the 96 Outlet wells
2. Lysis Mix to chambers 0+1+2 from Inlet #3
65°C, 5 min (300 sec)
10°C, 1 min (60 sec)
3. RT Mix to chambers 0+1+2+3 from Inlet #4
42°C, 120 min (7200 sec)
10°C, 1 min (60 sec)
4. 2nd Strand Mix to chambers 0+1+2+3+4 from Inlet #7
16°C, 120 min (7200 sec)
65°C, 20 min (1200 sec)
10°C, 1 min (60 sec)
5. IVT Mix to chambers 0+1+2+3+4+5 from Inlet #8
37°C, 12 hr (12× 3600 sec)
10°C, 1 min (60 sec)
6. Harvest

aRNA Clean up

- 1) Pool 2.5 µL of each of the harvested aRNA libraries (~240 µL total)
- 2) Add 1.8× volumes (430 µL) of RNAClean XP beads to the harvested aRNA. Mix well until the liquid appears homogeneous. Incubate for 10 min at room temperature.
- 3) Bind beads on a magnetic stand for 5 min.
- 4) Discard the supernatant.
- 5) Wash 3 times with 70% ethanol.
- 6) Air dry for 10 min or until completely dry.
- 7) Resuspend the beads with 16 µl RNase-free water. Incubate in room temperature for 2 min.
- 8) Bind beads on a magnetic stand for 5 min.
- 9) Transfer the eluted aRNA to a new Eppendorf, continue with ExoSAP treatment, as in regular CEL-Seq2 protocol.

EXO-SAP treatment (to remove primers):

- Add 6 µL enzyme
- Incubate at 15 minutes at 37 °C

RNA fragmentation:

- Mix the following on ice:

aRNA	22 μ L
Fragmentation buffer	5.5 μ L
- Incubate for 3 min. at 94°C.
- Immediately move to ice and add 2.75 μ L Fragmentation stop buffer.


aRNA cleanup:

- Prewarm RNAClean XP beads to room temperature.
- Vortex RNAClean XP beads until well dispersed, add to sample 55 μ L beads. (1.8 volumes)
- Incubate at room temperature for 10 min.
- Place on magnetic stand for at least 5 min, until liquid appears clear.
- Remove and discard 50 μ L of the supernatant.
- Add 200 μ L freshly prepared 70% EtOH.
- Incubate at least 30 seconds, then remove and discard supernatant without disturbing beads.
- Repeat wash two more times.
- Air dry beads for 15 min, or until completely dry.
- Resuspend with 7 μ L water. Pipette entire volume up and down ten times to mix thoroughly.
- Incubate at room temperature for 2 min.
- Place on magnetic stand for 5 min, until liquid appears clear.
- Transfer supernatant to new tube.

Stopping point: Samples can be kept at -80°C

Check aRNA amount and quality:

- Load 1 μ L onto Bioanalyzer RNA pico chip after heating an aliquot of the sample to 70° for 2 min.
- When starting the IVT with ~0.1 ng total RNA, the expected yield is 500-1000 pg/ μ L. Size distribution should peak at ~500 bp (See Bioanalyzer plot for example). If yield is higher, dilute a portion of aRNA to 1 ng/ μ L for the RT step below.



Library preparation:

RT reaction:

To 5 μL aRNA, add 1 μL randomhexRT primer and 0.5 μL 10 mM each dNTP.

Incubate 5 min at 65 $^{\circ}\text{C}$, quick chill on ice.

- Add 4 μL of the following mix at room temperature to each reaction:

First Strand buffer	2 μL
DTT 0.1M	1 μL
RNaseOUT	0.5 μL
SuperscriptII	0.5 μL

Incubate 10 min at 25 $^{\circ}\text{C}$.

Incubate 1 hr at 42 $^{\circ}\text{C}$ (in hybridization oven, or pre heated thermal cycler with lid at 50 $^{\circ}\text{C}$)

PCR amplification:

To each reverse transcription reaction add 38 μL of the following mix:

- Ultra Pure Water 11 μL
- PCR Master Mix 25 μL
- RNA PCR Primer (RP1, from Illumina kit) 2 μL

To each reaction add 2 μL of a uniquely indexed RNA PCR Primer (e.g. RPi1; if mixing multiple libraries for sequencing, choose balanced primers according to Illumina's pooling guide)

Amplify the tube in the thermal cycler using the following PCR cycling conditions:

- 30 seconds at 98°C
- 11 cycles of:
 - 10 seconds at 98°C
 - 30 seconds at 60°C
 - 30 seconds at 72°C
- 10 minutes at 72°C
- Hold at 4°C

Can go up to 15 cycles if necessary, if aRNA concentration was low.

Stopping point: samples can be kept at -20°C.

Bead Cleanup of PCR products – Repeat 1:

- Prewarm beads to room temperature.
- Vortex AMPure XP Beads until well dispersed, then add 50 µL to the 50 µL PCR reaction. Mix entire volume up ten times to mix thoroughly.
- Incubate at room temperature for 15 min.
- Place on magnetic stand for at least 5 min, until liquid appears clear.
- Remove and discard 95 µL of the supernatant.
- Add 200 µL freshly prepared 80% EtOH.
- Incubate at least 30 seconds, then remove and discard supernatant without disturbing beads.
- Add 200 µL freshly prepared 80% EtOH
- Incubate at least 30 seconds, then remove and discard supernatant without disturbing beads.
- Air dry beads for 15 min, or until completely dry.
- Resuspend with 25 µL water. Pipette entire volume up and down ten times to mix thoroughly.
- Incubate at room temperature for 2 min.
- Place on magnetic stand for 5 min, until liquid appears clear.
- Transfer 25 µL of supernatant to new tube.


Bead Cleanup of PCR products – Repeat 2:

Repeat as above, adding 25 µL beads and eluting in 10 µL water at the end, transferring 10 µL to a new tube.

Check library amount and quality:

Check concentration of DNA by Qubit, 1 µL should be enough to measure using the high sensitivity reagent; expected concentration is at least ~1ng/µL.

Run 1 µL of each sample on Bioanalyzer using a high sensitivity DNA chip to see size distribution. Expected peak at 200-400 bp (See Bioanalyzer plot for example).



Concentration to be loaded for sequencing should be calibrated by the sequencing facility. We are using 8 pM with HiSeq high throughput v.3 reagents, or 12 pM on MiSeq and HiSeq rapid mode. Paired end sequencing is performed, 15 bases for read 1, 7 for the illumine index (when needed) and 36 bases for read 2. Throughout the Illumina sequencing the libraries should be considered Small-RNA libraries. For example, in Illumina Rehyb. kits, some Illumina kits are not for Small-RNA.

Currently, CEL-Seq2 libraries are not compatible with Illumina HiSeq high throughput v.4 reagents.