


Name C1 CAGE
 Revision B
 Description C1 CAGE is a method for single-cell transcriptome analysis for molecular counting of RNA 5'-ends. Paired-end sequencing, random priming and unique molecular identifiers are used for single-molecule fragment assembly of mRNAs and long non-coding RNAs.
 Authors T.Kouno, S.Kato, M.Mendez, I.Abugessaisa, J.Shin and C.Plessy
 Institution RIKEN, JAPAN
 Lab Division of Genomics Technologies, Center for Life Sciences Technologies
 Special Instructions




Script Summary - Prime			
Runtime Estimates			
Barcode	Estimate		
1861x (5-10 um diameter cells)	0 hours, 11 minutes		
1862x (10-17 um diameter cells)	0 hours, 13 minutes		
1863x (17-25 um diameter cells)	0 hours, 12 minutes		
Script Summary - Cell Load			
Runtime Estimates			
Barcode	Estimate		
1861x (5-10 um diameter cells)	0 hours, 31 minutes		
1862x (10-17 um diameter cells)	1 hours, 0 minutes		
1863x (17-25 um diameter cells)	0 hours, 51 minutes		
Incubation Profile			
Script Step	Operation	Temperature (C)	Duration (s)
Stain	Incubation	25	600
Script Summary - Sample Prep			
Runtime Estimates			
Barcode	Estimate		
1861x (5-10 um diameter cells)	8 hours, 52 minutes		
1862x (10-17 um diameter cells)	8 hours, 52 minutes		
1863x (17-25 um diameter cells)	8 hours, 52 minutes		
Incubation Profile			
Script Step	Operation	Temperature (C)	Duration (s)
Lysis	Incubation	S1 72	180




		S2	4	600
		S3	25	60
RT	Incubation	S1	22	600
		Reverse Transcription	42	5400
		Enzyme inactivation	75	900
		S4	25	10
PCR	Hot Start	95C	95	60
PCR	PCR x30	Denaturation	95	15
		Annealling	65	30
		Extension	68	360
PCR	Extension	72C	72	600



Script Reagent Details - Prime



Reagent Loading			
Name	Volume (µl)	IFC Inlet	Notes
● Harvest Reagent	200 µl	A1	
● Harvest Reagent	200 µl	A2	
● Blocking Reagent	15 µl	C1	
● Blocking Reagent	15 µl	C2	
● Harvest Reagent	20 µl	P1	
● Harvest Reagent	20 µl	P2	
Reagent Mix Recipe - Prime			
Blocking Reagent			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Blocking RGT (1X)			
Harvest Reagent			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Harvest RGT (1X)			


Script Reagent Details - Cell Load


Inlet Reuse			
Name	IFC Inlet	Instructions	
● Cell Mix	C1	Aspirate inlet prior to loading reagents	
● Flowthrough Outlet	C2	Aspirate inlet prior to loading reagents (1862x, 1863x only)	
Reagent Loading			
Name	Volume (μl)	IFC Inlet	Notes
● LIVE/DEAD Staining Solution	7	1	
● Preloading Reagent	↻ 24	2	
● Cell Wash Buffer	7	5	
● Cell Mix	↻ 6	C1	
Reagent Mix Recipe - Cell Load			
Preloading Reagent			
Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Preloading RGT (1X)			
Cell Mix			
Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
Suspension RGT (2.5X)	40	1	1
Cells 250 / μL	60		
100 Total Prep Volume			
LIVE/DEAD Staining Solution			
Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
Cell Wash BUF (1X)	1250	0.9975	0.9975
Ethidium homodimer-1 (2 mM)	2.5	0.004	0.004
Calcein AM (4 mM)	0.625	0.002	0.002



1253.125 Total Prep Volume

Cell Wash Buffer			
Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
Cell Wash BUF (1X)			

Script Reagent Details - Sample Prep


Inlet Reuse				
Name	IFC Inlet	Instructions		
● Preloading Reagent	2	Aspirate inlet prior to loading reagents		
● Flowthrough Outlet	C2	Aspirate inlet prior to loading reagents (1862x, 1863x only)		
Reagent Loading				
Name	Volume (μl)	IFC Inlet	Notes	
● Preloading Reagent	24	2		
● C1 CAGE Lysis Mix with Diluted RNA Spikes	7	3		
● C1 CAGE RT Mix	8	4		
● C1 CAGE PCR Mix	24	7		
● C1 CAGE PCR Mix	24	8		
● Harvest Reagent	180 μl each	Harvest Inlets		
Reagent Mix Recipe - Sample Prep				
Preloading Reagent				
Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.	
C1 Preloading RGT (1X)				
C1 CAGE Diluted RNA spikes (Secondary)				
Special Instructions:				

Dilution of spike is performed by serial dilution.				
Like 1 uL ERCC RNA Spike-In Mix I to 9 uL C1 Loading RGT, and put 1 uL into 132.2 uL C1 Loading RGT.				
Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.	
C1 Loading RGT (20X)	1332	19.985		
ERCC RNA Spike-In Mix I	1			

1333 Total Prep Volume

C1 CAGE Lysis Mix with Diluted RNA Spikes

Special Instructions:

RT primer sequence: TCGTCGGCAGCGTCAGATGTGNNNNNN

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
Triton X-100 solution -10% in H2O (10%)	0.67	0.2975	0.1984
DNA suspension buffer	18.65		
RNase inhibitor (40 U/μl)	0.57	1.0124	0.6753
RT primer (SupN6i) (50 μM)	1.13	2.5089	1.6734
C1 CAGE Diluted RNA spikes	1.5		

22.52 Total Prep Volume

C1 CAGE RT Mix

Special Instructions:

TS oligo sequence: TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGNNNNNNNTATA(rG)(rG)(rG)

* r is RNA

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Loading reagent	1.05		
5x First strand buffer (5X)	9.8	1.6752	0.9565
0.1 M DTT (0.1M)	4.9	0.0168	0.0096
10mM dNTP Mix (10 mM)	3.08	1.053	0.6013
Betaine (5M)	3.08	0.5265	0.3006
RNase inhibitor (40 U/μl)	1.24	1.6957	0.9683
TS oligo (NC2lg_Nbi) (500 μM)	0.7	11.9658	6.8325
Superscript III (200 U/μl)	4.9	33.5043	19.1309
Nuclease-Free Water	0.5		

29.25 Total Prep Volume

C1 CAGE PCR Mix

Special Instructions:

PCR primer sequence: TCGTCGGCAGCGTCAGATGTG

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Loading reagent	3.75		
PCR-Grade water	54.75		
10x Advantage2 PCR buffer (not SA)	7.5		
50x dNTP Mix (10 mM)	3	0.4	0.3244
PCR primer (Supi) (12 μM)	3	0.48	0.3893
50x Advantage2 Polymerase Mix	3		

75 Total Prep Volume

Harvest Reagent



Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Harvest RGT (1X)			

Protocol Reagent Shopping List							
Reagent Name	Vendor	Part Number	Kit Part Number	Stock Concentration			
Ethidium homodimer-1	Life Technologies		L-3224	2 mM			
Calcein AM	Life Technologies		L-3224	4 mM			
Triton X-100 solution -10% in H2O	SIGMA	93443-100ML		10%			
DNA suspension buffer	TEKnova	T0221					
RNase inhibitor	TaKaRa	2313A	6110A	40 U/μl			
RT primer (SupN6i)	Invitrogen			50 μM			
ERCC RNA Spike-In Mix I	Ambion	4456740					
5x First strand buffer	Life Technologies	Y02321	18080044	5X			
0.1 M DTT	Invitrogen	Y00147		0.1M			
10mM dNTP Mix	TaKaRa		6110A	10 mM			
Betaine	Wako	027-10865		5M			
TS oligo (NC2lg_Nbi)	IDT			500 μM			
Superscript III	Life Technologies	56575	18080044	200 U/μl			
Nuclease-Free Water	Ambion	AM9937					
PCR-Grade water	TaKaRa	S1785	639206				
10x Advantage2 PCR buffer (not SA)	TaKaRa	S1799	639206				
50x dNTP Mix	TaKaRa	S1800	639206	10 mM			
PCR primer (Supi)	Invitrogen			12 μM			
50x Advantage2 Polymerase Mix	TaKaRa	S1798	639206				
Fluidigm Reagent Kits							
Reagent Name	Part Number	Stock Concentration	PN 100-8920	PN 100-6201	PN 100-5319	PN 100-7357	PN 100-8921
C1 Blocking RGT	100-5316	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Harvest RGT	100-6248	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	
C1 Preloading RGT	100-5311	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Suspension RGT	100-5315	2.5X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Cell Wash BUF	100-5314	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
C1 Loading RGT	100-5170	20X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
C1 Loading reagent	100-5170		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>