


Name CEL Seq 2  
 Revision A  
 Description CELseq2 is a 3'-end counting mRNAseq method that uses in vitro transcription for initial amplification. CELseq2 incorporates numerous improvements over the original method including increased sensitivity and incorporation of unique molecular identifiers.  
 Authors Tamar Hashimshony and Itai Yanai  
 Institution Technion - Israel Institute of Technology  
 Lab Itai Yanai  
 Special Instructions T. Hashimshony, T. N. Senderovich, G. Avital, A. Klochendler, Y. de Leeuw, L. Anavy, D. Gennery, S. Li, K. J. Livak, O. Rozenblatt-Rosen, Y Dor, A. Regev, and I. Yanai. CEL-Seq2: sensitive highly-multiplexed single-cell RNA-Seq. Genome Biol 17:77 (2016). Complete protocol including sequences of all primers and post-C1 processing procedure are in the document entitled "CELSeq2 Protocol for Script Hub"



### Script Summary - Prime

#### Runtime Estimates

Barcode	Estimate
1861x (5-10 um diameter cells)	0 hours, 11 minutes
1862x (10-17 um diameter cells)	0 hours, 13 minutes
1863x (17-25 um diameter cells)	0 hours, 12 minutes

### Script Summary - Capture

#### Runtime Estimates

Barcode	Estimate
1861x (5-10 um diameter cells)	0 hours, 15 minutes
1862x (10-17 um diameter cells)	0 hours, 34 minutes
1863x (17-25 um diameter cells)	0 hours, 27 minutes

### Script Summary - CELseq

#### Runtime Estimates

Barcode	Estimate
1861x (5-10 um diameter cells)	19 hours, 55 minutes
1862x (10-17 um diameter cells)	19 hours, 55 minutes
1863x (17-25 um diameter cells)	19 hours, 55 minutes

#### Incubation Profile


Script Step	Operation	Temperature (C)	Duration (s)
Backload RT primers	Incubation	25	300
Lysis	Lysis	65	300



		Cooling	10	60
RT	Reverse Transcriptase Reaction	Incubation	42	7200
		Cooling	10	60
2nd Strand	2nd Strand Synthesis	Incubation	16	7200
		Enzyme Inactivation	65	1200
		Cooling	10	60
IVT	T7 Transcription x12	Incubation	37	3600
IVT	Cooling	Cooling	10	60




**Script Reagent Details - Prime**



Reagent Loading			
Name	Volume (µl)	IFC Inlet	Notes
● Harvest Reagent	200 µl	A1	
● Harvest Reagent	200 µl	A2	
● Blocking Reagent	15 µl	C1	
● Blocking Reagent	15 µl	C2	
● Harvest Reagent	20 µl	P1	
● Harvest Reagent	20 µl	P2	
Reagent Mix Recipe - Prime			
Blocking Reagent			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Blocking RGT (1X)			
Harvest Reagent			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Harvest RGT (1X)			



## Script Reagent Details - Capture




Inlet Reuse			
Name	IFC Inlet	Instructions	
● Cell Mix	C1	Aspirate inlet prior to loading reagents	
● Flowthrough Outlet	C2	Aspirate inlet prior to loading reagents (1862x, 1863x only)	
Reagent Loading			
Name	Volume (µl)	IFC Inlet	Notes
● Preloading Reagent	20	2	
● Cell Wash Buffer	7	5	
● Cell Mix	6	C1	
Reagent Mix Recipe - Capture			
Preloading Reagent			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Preloading RGT (1X)			
Cell Mix			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Suspension RGT (2.5X)	40	1	1
Cells 66-330 / µL	60		
100 Total Prep Volume			
Cell Wash Buffer			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Cell Wash BUF (1X)			




**Script Reagent Details - CELseq**

Sample Prep Description



Inlet Reuse			
Name	IFC Inlet	Instructions	
● Bus Wash	2	Aspirate inlet prior to loading reagents	
● Flowthrough Outlet	C2	Aspirate inlet prior to loading reagents (1862x, 1863x only)	
Reagent Loading			
Name	Volume (µl)	IFC Inlet	Notes
● Bus Wash	20	2	
● Lysis mix	7	3	
● RT mix	7	4	
● 2nd Strand mix	24	7	
● IVT mix	24	8	
● Harvest Reagent	180 µl each	Harvest Inlets	
● CEL-Seq2 primers	5 each	Harvest Outlets	
Reagent Mix Recipe - CELseq			
Bus Wash			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
Not Defined			
5X Lysis Buffer (Secondary: 5X)			
Special Instructions:			
-----			
Can be made up as a 100 ml stock by using ml in the recipe rather than µl. Store at room temperature.			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
NP-40 (10%)	25	2.5	
Tris-HCl, pH 8.4 (1000 mM)	25	250	



EDTA (500 mM)	1	5
Water	49	

100 Total Prep Volume

**Lysis mix**

Special Instructions:

-----

1 ml is made so there is enough Lysis Mix to dilute the CEL-Seq2 primers (see previous step).

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
RNaseOUT (40 U/μl)	4	0.16	0.064
C1 Loading Reagent (20X)	70	1.4	0.56
Water	610		
1:10,000 dilution of ERCC RNA Spike-In Mix (0.01%)	8	0.0001	0
5X Lysis Buffer (5X)	280	1.4	0.56
DTT (dithiothreitol) (1000)	28	28	11.2

1000 Total Prep Volume

**RT mix**

Comments:

-----

First Strand Buffer: included with SuperScript II

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
First Strand Buffer (5X)	20	3.4965	1
dNTP Mix (each) (25 mM)	2	1.7483	0.5
RNaseOUT (40 U/μl)	2.6	3.6364	1.04
SuperScript II Reverse Transcriptase (200 U/μl)	2.6	18.1818	5.2
C1 Loading Reagent (20X)	1.4	0.979	0.28

28.6 Total Prep Volume

**IVT mix**

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
T7 Reaction Buffer (10X)	6.6	2.2297	0.9989
T7 Enzyme Mix (10X)	6.6	2.2297	0.9989
C1 Loading Reagent (20X)	1.5	1.0135	0.4541
rNTP Mix (each) (25 mM)	14.9	12.5845	5.6378

29.6 Total Prep Volume

**Harvest Reagent**

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Harvest Reagent (1X)	180	1	

180 Total Prep Volume

**CEL-Seq2 primers**

Special Instructions:

9/12/2016 11:04:45 PM



-----

Each CEL-Seq2 primer is one of 96, each with its own sample barcode. Sequences are in the document "CEL-Seq2 Protocol for Script Hub." These primer stocks should be stored in a 96-well plate at -20 degrees. The process is to dispense 8 µl Lysis Mix (see Step: Lysis) to each well of a 96-well plate. Then 2 µl primer is added to each well from the primer stock plate using an 8-channel pipettor. The 96-well plate with Lysis Mix plus primers is then gently vortexed and briefly centrifuged.

Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
Lysis Mix (1X)	8	0.8	0.5336
CEL-Seq2 Primer (25 µg/ml)	2	5	3.335

10 Total Prep Volume

<b>2nd Strand mix</b>			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
Second Strand Buffer (5X)	24.7	1.235	1.0016
dNTP Mix (each) (25 mM)	1	0.25	0.2028
E coli DNA Polymerase I (10 U/µl)	3.3	0.33	0.2676
RNase H (2 U/µl)	1	0.02	0.0162
C1 Loading Reagent (20X)	5	1	0.811
Water	64		
E coli Ligase (10 U/µl)	1	0.1	0.0811

100 Total Prep Volume

**Protocol Reagent Shopping List**

Reagent Name	Vendor	Part Number	Kit Part Number	Stock Concentration
RNaseOUT	Thermo Fisher	10777019		40 U/μl
1:10,000 dilution of ERCC RNA Spike-In Mix	Ambion	4456740		0.01%
NP-40	Thermo Fisher	28324		10%
Tris-HCl, pH 8.4	Teknova	T1084		1000 mM
EDTA	Teknova	E0306		500 mM
DTT (dithiothreitol)	Teknova	D9750		1000
dNTP Mix (each)	Thermo Fisher	R1121		25 mM
SuperScript II Reverse Transcriptase	Thermo Fisher	18064014		200 U/μl
T7 Reaction Buffer	Thermo Fisher		AM1334	10X
T7 Enzyme Mix	Thermo Fisher		AM1334	10X
rNTP Mix (each)	New England BioLabs	N0466S		25 mM
Second Strand Buffer	Thermo Fisher	10812014		5X
E coli DNA Polymerase I	Thermo Fisher	18010025		10 U/μl
RNase H	Thermo Fisher	18021071		2 U/μl

**Fluidigm Reagent Kits**

Reagent Name	Part Number	Stock Concentration	PN 100-8920	PN 100-6201	PN 100-5319	PN 100-7357	PN 100-8921
C1 Blocking RGT	100-5316	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Harvest RGT	100-6248	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
C1 Preloading RGT	100-5311	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Suspension RGT	100-5315	2.5X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Cell Wash BUF	100-5314	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
C1 Loading Reagent	100-5170	20X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
C1 Harvest Reagent	100-7081	1X					