


Name	LPS Simulation
Revision	B
Description	script used in Shalek et al., 2014 Stimulation of cells with LPS
Authors	Shalek AK, Satija R, Shuga J, Trombetta JJ, Gennert D, et al.
Institution	Broad Institute, Harvard University, MIT, Hebrew University, Fluidigm Corporation
Lab	Regev, Park and May
Special Instructions	



Script Summary - Prime

Runtime Estimates

Barcode	Estimate
1861x (5-10 um diameter cells)	0 hours, 11 minutes
1862x (10-17 um diameter cells)	0 hours, 13 minutes
1863x (17-25 um diameter cells)	0 hours, 12 minutes
1771x (5-10 um diameter cells)	0 hours, 11 minutes
1772x (10-17 um diameter cells)	0 hours, 13 minutes
1773x (17-25 um diameter cells)	0 hours, 12 minutes

Script Summary - Cell Load and Stimulation

Runtime Estimates

Barcode	Estimate
1861x (5-10 um diameter cells)	0 hours, 25 minutes
1862x (10-17 um diameter cells)	0 hours, 44 minutes
1863x (17-25 um diameter cells)	0 hours, 37 minutes
1771x (5-10 um diameter cells)	0 hours, 25 minutes
1772x (10-17 um diameter cells)	0 hours, 44 minutes
1773x (17-25 um diameter cells)	0 hours, 37 minutes


Incubation Profile

Script Step	Operation	Temperature (C)	Duration (s)
Stimulation	Incubation	37	600


Script Summary - Sample Prep

Runtime Estimates


Barcode	Estimate
1861x (5-10 um diameter cells)	8 hours, 6 minutes



1862x (10-17 um diameter cells)	8 hours, 16 minutes			
1863x (17-25 um diameter cells)	8 hours, 14 minutes			
1771x (5-10 um diameter cells)	8 hours, 6 minutes			
1772x (10-17 um diameter cells)	8 hours, 16 minutes			
1773x (17-25 um diameter cells)	8 hours, 14 minutes			
Incubation Profile				
Script Step	Operation		Temperature (C)	Duration (s)
Lysis	Incubation	S1	72	180
		S2	4	600
		S3	25	60
RT	Incubation	S1	42	5400
		S2	70	600
PCR	Hot Start	95C	95	60
PCR	PCR x5	Denaturation	95	20
		Annealling	58	240
		Extension	68	360
PCR	PCR x9	Denaturation	95	20
		Annealling	64	30
		Extension	68	360
PCR	PCR x7	Denaturation	95	30
		Annealling	64	30
		Extension	68	420
PCR	Extension	72C	72	600



Script Reagent Details - Prime



Reagent Loading

Name	Volume (μl)	IFC Inlet	Notes
● Harvest Reagent	200 μl	A1	
● Harvest Reagent	200 μl	A2	
● Blocking Reagent	15 μl	C1	
● Blocking Reagent	15 μl	C2	
● Harvest Reagent	20 μl	P1	
● Harvest Reagent	20 μl	P2	


Reagent Mix Recipe - Prime

Blocking Reagent

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Blocking RGT (1X)			


Harvest Reagent

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Harvest RGT (1X)			



Script Reagent Details - Cell Load and Stimulation

Cell load and LPS stimulation



Inlet Reuse

Name	IFC Inlet	Instructions
● Cell Mix	C1	Aspirate inlet prior to loading reagents
● Flowthrough Outlet	C2	Aspirate inlet prior to loading reagents (1862x, 1863x only)

Reagent Loading

Name	Volume (μl)	IFC Inlet	Notes
● Medium with LPS	7	1	
● Preloading Reagent	20	2	
● Cell Mix	6	C1	

Reagent Mix Recipe - Cell Load and Stimulation

Preloading Reagent

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Preloading RGT (1X)			

Cell Mix

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Suspension RGT (2.5X)	40	1	1
Cells 66 -330 / μL	60		

100 Total Prep Volume


Medium with LPS

Special Instructions:


Incubate in tissue culture incubator with environment control for desired time

Note: any activators or inhibitors can also be used.


Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
Cell culture medium (1X)	9990	0.999	0.999



Ultrapure LPS-EK (100 mg/ml)	10	0.1	0.1
10000 Total Prep Volume			



Script Reagent Details - Sample Prep



Inlet Reuse

Name	IFC Inlet	Instructions
● Preloading Reagent	2	Aspirate inlet prior to loading reagents
● Flowthrough Outlet	C2	Aspirate inlet prior to loading reagents (1862x, 1863x only)

Reagent Loading

Name	Volume (μl)	IFC Inlet	Notes
● Preloading Reagent	20	2	
● RNA Seq Lysis Mix with Diluted RNA Spikes	7	3	
● RNA Seq RT Mix	8	4	
● Cell Wash Buffer	7	5	
● RNA Seq PCR Mix	24	7	
● RNA Seq PCR Mix	24	8	
● Harvest Reagent	180 μl each	Harvest Inlets	

Reagent Mix Recipe - Sample Prep

Preloading Reagent

Reagent (Stock Concentration)	Mix Prep (μl)	Prep Conc.	Chamber Conc.
C1 Preloading RGT (1X)			

ArrayControl RNA Spikes (Secondary: 1X)


Special Instructions:

After the ArrayControl RNA Spikes have thawed, remove spikes 1, 4, and 7 from the box. Pipette the following in three tubes:

Tube A - 13.5 μl THE RNA Storage Solution, RNA Spikes #7 - 1.5 μl

Tube B - 12.5 μl THE RNA Storage Solution, RNA Spikes #4 - 1.5 μl

Tube C - 148.5 μl THE RNA Storage Solution, RNA Spikes #1 - 1.5 μl



Vortex briefly tube A and spin to collect contents. Pipette 1.5 µl from tube A into tube B. Discard tube A.

Vortex briefly tube B and spin to collect contents. Pipette 1.5 µl from tube B into tube C. Discard tube B.
Vortex briefly tube C and spin to collect contents.

Tube C is the concentrated RNA Standard that may be aliquoted and frozen for future use.
Aliquot in tubes containing 1 µl volumes and store at -80 °C until use. One tube is necessary for each C1 chip run.

Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
THE RNA Storage Solution	174		
ArrayControl RNA Spikes #1 (100 ng/µl)	1.5	0.8403	
ArrayControl RNA Spikes #4 (100 ng/µl)	1.5	0.8403	
ArrayControl RNA Spikes #7 (100 ng/µl)	1.5	0.8403	

178.5 Total Prep Volume

RNA Seq Diluted RNA spikes (Secondary: 20X)

Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Loading RGT (20X)	99	19.8	
ArrayControl RNA Spikes (1X)	1	0.01	

100 Total Prep Volume

RNA Seq Lysis Mix with Diluted RNA Spikes

Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
RNA Seq Diluted RNA spikes (20X)	1	1	0.667
RNase Inhibitor (40 U/µl)	0.5	1	0.667
3' SMART CDS Primer IIA (12 µM)	7	4.2	2.8014
Clontech Dilution Buffer (1X)	11.5	0.575	0.3835

20 Total Prep Volume


RNA Seq RT Mix

Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Loading RGT (20X)	1.2	0.75	0.4283
5X First-Strand Buffer (Rnase-free) (5X)	11.2	1.75	0.9992
Dithiothreitol (100 mM)	1.4	4.375	2.4981
dNTP Mix (10 mM)	5.6	1.75	0.9992
SMARTer IIA Oligonucleotide (12 µM)	5.6	2.1	1.1991
RNase Inhibitor (40 U/µl)	1.4	1.75	0.9992
SMARTScribe Reverse Transcriptase (100 U/µl)	5.6	17.5	9.9925

32 Total Prep Volume

RNA Seq PCR Mix

Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
-------------------------------	---------------	------------	---------------




C1 Loading RGT (20X)	4.5	1	0.811
PCR-Grade Water	63.5		
10X Advantage 2 PCR Buffer (10X)	10	1.1111	0.9011
50X dNTP Mix (10 mM)	4	0.4444	0.3604
IS PCR primer (12 µM)	4	0.5333	0.4325
50X Advantage 2 Polymerase Mix (50X)	4	2.2222	1.8022

90 Total Prep Volume

Harvest Reagent			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Harvest RGT (1X)			
Cell Wash Buffer			
Special Instructions:			

Wash step: Wash cells after incubation and before mRNA preparation			
Reagent (Stock Concentration)	Mix Prep (µl)	Prep Conc.	Chamber Conc.
C1 Cell Wash BUF (1X)	20	1	1

20 Total Prep Volume

**Protocol Reagent Shopping List**

Reagent Name	Vendor	Part Number	Kit Part Number	Stock Concentration
Ultrapure LPS-EK	Invivogen	TLRL-PEKLPS		100 mg/ml
THE RNA Storage Solution	Life Technologies	AM7000		
ArrayControl RNA Spikes #1	Life Technologies		AM1780	100 ng/μl
ArrayControl RNA Spikes #4	Life Technologies		AM1780	100 ng/μl
ArrayControl RNA Spikes #7	Life Technologies		AM1780	100 ng/μl
RNase Inhibitor	Clontech		634833	40 U/μl
3' SMART CDS Primer IIA	Clontech		634833	12 μM
Clontech Dilution Buffer	Clontech			1X
5X First-Strand Buffer (Rnase-free)	Clontech		634833	5X
Dithiothreitol	Clontech		634833	100 mM
dNTP Mix	Clontech		634833	10 mM
SMARTer IIA Oligonucleotide	Clontech		634833	12 μM
SMARTScribe Reverse Transcriptase	Clontech		634833	100 U/μl
PCR-Grade Water	Clontech		639207	
10X Advantage 2 PCR Buffer	Clontech		639207	10X
50X dNTP Mix	Clontech		639207	10 mM
IS PCR primer	Clontech		634833	12 μM
50X Advantage 2 Polymerase Mix	Clontech		639207	50X

Fluidigm Reagent Kits

Reagent Name	Part Number	Stock Concentration	PN 100-8920	PN 100-6201	PN 100-5319	PN 100-7357	PN 100-8921
C1 Blocking RGT	100-5316	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Harvest RGT	100-6248	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	
C1 Preloading RGT	100-5311	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Suspension RGT	100-5315	2.5X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
C1 Loading RGT	100-5170	20X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>			<input checked="" type="checkbox"/>
C1 Cell Wash BUF	100-5314	1X	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		