




Name sc-GEM
 Revision A
 Description sc-GEM
 Authors Lih Feng Cheow
 Institution Institute of Molecular and Cell Biology
 Lab Microfluidics Systems Biology Lab
 Special Instructions

Cell Load

Sample Prep



Script Summary - Prime

Runtime Estimates

| Barcode | Estimate |
|---------------------------------|---------------------|
| 1861x (5-10 um diameter cells) | 0 hours, 11 minutes |
| 1862x (10-17 um diameter cells) | 0 hours, 12 minutes |
| 1863x (17-25 um diameter cells) | 0 hours, 12 minutes |

Script Summary - Cell Load

Runtime Estimates

| Barcode | Estimate |
|---------------------------------|---------------------|
| 1861x (5-10 um diameter cells) | 0 hours, 15 minutes |
| 1862x (10-17 um diameter cells) | 0 hours, 30 minutes |
| 1863x (17-25 um diameter cells) | 0 hours, 27 minutes |


Script Summary - Sample Prep

Runtime Estimates


| Barcode | Estimate |
|---------------------------------|---------------------|
| 1861x (5-10 um diameter cells) | 11 hours, 8 minutes |
| 1862x (10-17 um diameter cells) | 11 hours, 8 minutes |
| 1863x (17-25 um diameter cells) | 11 hours, 8 minutes |

Incubation Profile


| Script Step | Operation | | Temperature (C) | Duration (s) |
|-------------|-----------|-----------|-----------------|--------------|
| Lysis | Lysis | Manual | 25 | 300 |
| | | lysis | 70 | 600 |
| | | New Slice | 25 | 30 |
| RT | RT | S1 | 25 | 600 |
| | | S2 | 42 | 3600 |
| | | New Slice | 70 | 900 |



| | | | | |
|----------|------------------------|--------------|----|------|
| | | New Slice | 25 | 30 |
| Protease | Protease | Protease | 25 | 10 |
| | | New Slice | 50 | 5400 |
| | | New Slice | 70 | 1800 |
| | | New Slice | 25 | 30 |
| RD | RD | S1 | 25 | 10 |
| | | New Slice | 37 | 7200 |
| | | New Slice | 80 | 1200 |
| | | New Slice | 25 | 30 |
| PreAmp | Pre-Amp | Manual | 95 | 600 |
| PreAmp | PCR Cycle (3 Step) x22 | Denaturation | 95 | 15 |
| | | Annealling | 60 | 240 |




Script Reagent Details - Prime




| Reagent Loading | | | |
|-------------------------------|---------------|------------|---------------|
| Name | Volume (µl) | IFC Inlet | Notes |
| ● Harvest Reagent | 200 µl | A1 | |
| ● Harvest Reagent | 200 µl | A2 | |
| ● Blocking Reagent | 15 µl | C1 | |
| ● Blocking Reagent | 15 µl | C2 | |
| ● Harvest Reagent | 20 µl | P1 | |
| ● Harvest Reagent | 20 µl | P2 | |
| Reagent Mix Recipe - Prime | | | |
| Blocking Reagent | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Not Defined | | | |
| Harvest Reagent | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Not Defined | | | |




Script Reagent Details - Cell Load




| Inlet Reuse | | | |
|--------------------------------|---------------|--|---------------|
| Name | IFC Inlet | Instructions | |
| ● Cell Mix | C1 | Aspirate inlet prior to loading reagents | |
| ● Flowthrough Outlet | C2 | Aspirate inlet prior to loading reagents (1862x, 1863x only) | |
| Reagent Loading | | | |
| Name | Volume (µl) | IFC Inlet | Notes |
| ● Preloading Reagent | 20 | 2 | |
| ● Cell Wash Buffer | 7 | 5 | |
| ● Cell Mix | 5 | C1 | |
| Reagent Mix Recipe - Cell Load | | | |
| Preloading Reagent | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Not Defined | | | |
| Cell Mix | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Cells 66-330 / µL | | | |
| Cell Wash Buffer | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Cell Wash BUF (1X) | | | |



Script Reagent Details - Sample Prep



| Inlet Reuse | | | |
|----------------------------------|---------------|--|---------------|
| Name | IFC Inlet | Instructions | |
| ● Preloading Reagent | 2 | Aspirate inlet prior to loading reagents | |
| ● Flowthrough Outlet | C2 | Aspirate inlet prior to loading reagents (1862x, 1863x only) | |
| Reagent Loading | | | |
| Name | Volume (µl) | IFC Inlet | Notes |
| ● Preloading Reagent | 20 | 2 | |
| ● Lysis Final Mix | 7 | 3 | |
| ● RT | 7 | 4 | |
| ● Protease | 7 | 6 | |
| ● RD | 24 | 7 | |
| ● Pre Amp | 24 | 8 | |
| ● Harvest Reagent | 180 µl each | Harvest Inlets | |
| Reagent Mix Recipe - Sample Prep | | | |
| Preloading Reagent | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Not Defined | | | |
| Lysis Final Mix | | | |
| Reagent (Stock Concentration) | Mix Prep (µl) | Prep Conc. | Chamber Conc. |
| Single-Cell Lysis Solution (2) | 12.75 | 1.4167 | 0.9449 |
| C1 Lysis Plus RGT | 4.35 | | |
| PCR Water | 0.9 | | |
| 18 Total Prep Volume | | | |
| RT | | | |



| Reagent (Stock Concentration) | Mix Prep (μl) | Prep Conc. | Chamber Conc. |
|-------------------------------|---------------|------------|---------------|
| Stop solution (15) | 2.91 | 2.425 | 0.97 |
| Single-Cell VILO RT Mix (5) | 8.76 | 2.4333 | 0.9733 |
| C1 Loading RGT (20X) | 0.9 | 1 | 0.4 |
| Single-Cell PreAmp Mix (5) | 5.43 | 1.5083 | 0.6033 |

18 Total Prep Volume

| Protease | | | |
|-------------------------------|---------------|------------|---------------|
| Reagent (Stock Concentration) | Mix Prep (μl) | Prep Conc. | Chamber Conc. |
| C1 Loading RGT (20X) | 2 | 1 | 0.286 |
| PCR Water | 18 | | |
| QIAGEN Protease (5 mg/ml) | 20 | 2.5 | 0.715 |

40 Total Prep Volume

| Harvest Reagent | | | |
|-------------------------------|---------------|------------|---------------|
| Reagent (Stock Concentration) | Mix Prep (μl) | Prep Conc. | Chamber Conc. |
| C1 Harvest RGT (1X) | | | |

RD

| Reagent (Stock Concentration) | Mix Prep (μl) | Prep Conc. | Chamber Conc. |
|-------------------------------|---------------|------------|---------------|
| NEB Cutsmart Buffer (10X) | 4 | 1.1111 | 0.9011 |
| PCR Water | 28 | | |
| HpaII (50 Units/ul) | 4 | 5.5556 | 4.5056 |

36 Total Prep Volume


Pre Amp

Special Instructions:

To create the preamplification primer mix, all preamplification primers to be used should be combined and diluted in nuclease free water to a final concentration of 250 nM each primer.

| Reagent (Stock Concentration) | Mix Prep (μl) | Prep Conc. | Chamber Conc. |
|---|---------------|------------|---------------|
| C1 PreAmp Dilution RGT | 16 | | |
| Single-Cell PreAmp Mix (5) | 8 | 1 | 0.448 |
| Preamplification Primer mix (250 nM each primer) (250 nM) | 16 | 100 | 44.8 |

40 Total Prep Volume



Protocol Reagent Shopping List

| Reagent Name | Vendor | Part Number | Kit Part Number | Stock Concentration |
|----------------------------|---------------------|-------------|-----------------|---------------------|
| Single-Cell Lysis Solution | Life Technologies | | 4458237 | 2 |
| Stop solution | Life Technologies | | 4458237 | 15 |
| Single-Cell VILO RT Mix | Life Technologies | | 4458237 | 5 |
| Single-Cell PreAmp Mix | Life Technologies | | 4458237 | 5 |
| QIAGEN Protease | QIAGEN | 19155 | | 5 mg/ml |
| NEB Cutsmart Buffer | New England BioLabs | B7204S | | 10X |
| HpaII | New England BioLabs | R0171M | | 50 Units/ul |

Fluidigm Reagent Kits

| Reagent Name | Part Number | Stock Concentration | PN 100-5319 | PN 100-8920 | PN 100-6201 | PN 100-8921 | PN 100-7357 |
|------------------------|-------------|---------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| Cell Wash BUF | 100-5314 | 1X | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | | |
| C1 Lysis Plus RGT | 100-5313 | | <input checked="" type="checkbox"/> | | | | |
| PCR Water | 100-5941 | | | | | | <input checked="" type="checkbox"/> |
| C1 Loading RGT | 100-5170 | 20X | | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | <input checked="" type="checkbox"/> | |
| C1 Harvest RGT | 100-5310 | 1X | <input checked="" type="checkbox"/> | | | | |
| C1 PreAmp Dilution RGT | 100-5318 | | <input checked="" type="checkbox"/> | | | | |