

C1 MiniMetagenomics Qiagen Repli-g sc Protocol

Reagents and materials to UV treat (for 30min)
DLB2 (made from DLB+IM DTT)
H2O (high quality, PCR water)
PBS
1% NaCl
Blocking Solution for C1
Cell wash buffer (DNA seq)
Tris-Cl Ph 8.0 (10mM)
Repli-g stop solution
Preload Reagent for C1
1.5 mL Eppendorf tubes (about 6)
C1 chip

Other Reagents Not UV Treated
Harvest reagent for C1
10% tween
Lysozyme Stock
Repli-g Single Cell Polymerase
Repli-g Single Cell Reaction Buffer

Cell Sample Preparation:

Centrifuged quickly for the dirt to settle in the tube (typically using 2 mL tube)
Transfer clear sample into a clean tube
Spin down the sample at 5000g for 7 min
Resuspend in 1% UV treated NaCl
Dilute until ~2 cells per nL (2×10^6 cells / mL).

Prime the Chip:

10 uL of C1 blocking reagent in cell input and output ports
200 uL harvest reagent in top and bottom rubber sealed chambers
20 ul in 4 wells (2+2) on the left and right side of the chip
20 uL harvest reagent in 4*9 wells
Strip film from bottom of chip
Load chip
Run Mini-Metagenomics Prime Script

Load Cells:

Remove C1 blocking reagent
Put 10-15 uL of the properly diluted cell mixture only in the left cell load port
20 uL C1 pre-load reagent in well #2
Load chip
Run Mini-Metagenomics Random Capture Script

Making Master Mixes:

Make 10% tween in clean water

Make DLB2 solution

DLB 33 uL (44 ul)

DTT 1M 3 uL (4 ul)

UV treat for 30 min

Make Lysozyme Stock (Once):

Lysozyme Stock at 36,000 U/uL 2.2 uL

UV treated Tris-Cl (10mM) 197.9 uL

Final Lysozyme Concentration 400 U/uL

A	Make Lysozyme Enzyme Master Mix	
	400 U/uL Lysozyme Stock	10 uL
	UV treated Tris-Cl (10mM)	25 uL
	10% Tween	5 uL
	Final Lysozyme Concentration	100 U/uL
B	Make Lysis solution	
	DLB2 solution	19 uL
	10% Tween	1 uL
	Total :	20 uL
C	Make Stop Mix	
	Repli-g stop solution	19 uL
	10% Tween	1 uL
	Total:	20 uL
D	Make Salt Addition Solution with Master Mix	
	DLB2	2.25 uL
	Repli-g stop solution	2.25 uL
	Repli-g Single Cell Reaction Buffer	13 uL
	UV treated water	16.25 uL
	Total:	33.75 uL
E	Make MDA Master Mix	
	Repli-g Single Cell Reaction Buffer	30.75 uL
	Repli-g single cell polymerase	3 uL
	Total:	33.75 uL

Running MDA:

Add 7 uL lysozyme enzyme mix (A) to well #3

Add 7 uL lysis mix (B) to well #4

Add 9 uL Stop mix (C) to well #6

Add 24 uL (12 ul+12 ul) Make Salt Addition Solution with Master Mix (D) to well #7

Add 24 uL (12 ul+12 ul) Make MDA Master Mix (E) into well #8

Top up preload reagent to 20 uL in well #2

Add 180 uL harvest reagent to the 4 troughs

Load chip and run the Mini-Metagenomics Qiagen Lysis MDA script